



An ISO 9001 & 14001 Co.

MEDICAL & HOSPITAL EQUIPMENT



MODEL NO.:- RAJAS-CT2-16

(Venipuncture and injection arm (infusion, blood transfusion))

Specification:-

Venipuncture and injection arm (infusion, blood transfusion)

- Arm exercise including all venous access for intravenous treatment and incision of veins also sites for intramuscular injection and subcutaneous.
- With an extensive 8-branch vascular system or more allowing students to practice perform venipuncture at all sites adequate first and second rank, including including the suspension of infusions and the introduction of peripheral catheters.
- Supply of all veins with artificial blood from a single external fluid pocket.
- The back of the hand must be equipped with veins metacarpals and digital veins which can receive an injection.
- With sites prepared on the arm for injection intracutaneous
- Intramuscular injection into the M. deltoideus must be realistic in this area thanks to a skin soft and practically natural and by the points virtually natural bone orientation.
- Subcutaneous injections with distilled water must cause skin papules characteristics at the intended arm sites.

RAJAS ENTERPRISES (INDIA)

#59, Phase-1, Sector-1, HSIIDC, Industrial Estate, SAHA, Ambala Cantt-133104 (India)

Email:rajasmmedical@gmail.com, rajasglass@gmail.com

www.rajasenterprises.com