



An ISO 9001 & 14001 Co.

MEDICAL & HOSPITAL EQUIPMENT



MODEL NO.:- RAJAS-CT5-27
(Resuscitation table)

Specification:-

Resuscitation table

it must be equipped with radiant heating with a temperature control system: A heat source, often infrared, helps maintain the newborn's body temperature at an optimal level, often around 37°C. The bed must be reclining with folding side panels and equipped with an examination lamp and an LCD screen to display monitoring data, such as the baby's temperature, heart rate, etc. It must also be equipped with a phototherapy unit for the treatment of jaundice, a suction system, a drawer for storing equipment and a holder for infusions.

Delivery will be made to the following sites 11 by DRS DRS south-central (Manga) DRS central-plateau (Ziniaré), DRS north-central (KAYA), DRS EST (Fada) DRS CASCADES (Banfora) DRS north Ouahigouya), DRS CENTRE-EST (Tenkodogo), DRS Boucle du Mouhoun (Dédougou) DRS Centre Ouagadougou, DRS Hauts-Bassins (Bobo) DRS Sahel (DORI) South-west (Gaoua)

RAJAS ENTERPRISES (INDIA)

#59, Phase-1, Sector-1, HSIIDC, Industrial Estate, SAHA, Ambala Cantt-133104 (India)

Email: rajasmedical@gmail.com, rajasglass@gmail.com

www.rajasenterprises.com